PRODUCT-SPECIFIC INFORMATION ON THIS PAGE ONLY

Disposable Chemical Resistant Examination and Protective Gloves Chemstar® by Granberg. Nitrile, non-sterile, powder-free, orange.

AQL 0.65

EN ISO 21420:2020
ASTM D6978-05

Table with columns: Available sizes, S, M, L, XL, 2XL and rows: 6/7, 7/8, 8/9, 9/10, 10/11

Table with columns: EN ISO 374-1:2016+A1:2018 (Type B), Permeation Performance Level, Measured Breakthrough Time (minutes), EN ISO 374-4:2019 Mean Degradation, %

Chemotherapy Drugs tested in accordance with ASTM D6978-05.

Table with columns: Chemotherapy Drug in accordance with ASTM D6978-05, Minimum breakthrough detection time in minutes

Latex free: yes.

This product is Category III Personal Protective Equipment as per Regulation (EU) 2016/425 and complies with standards: EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016.

Notified Body responsible for EU Type Examination (Module B) and Quality Assurance of the Production Process (Module D): BSI Group The Netherlands B.V. (NB No. 2797), Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands.

This product is classified as Class I Medical Device according to Annex VIII of the Regulation (EU) 2017/745 and complies with standards: EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 15223-1:2021.

EU Declaration of Conformity: www.granberg.no/search

Head office: GRANBERG AS, Bjoavegen 1442, 5584 Bjoa, Norway. Phone: +47 53 77 53 00 E-mail: post@granberg.no

Swedish office: GRANBERG SVERIGE AB, Schubergsvägen 20, 311 74 Falkenberg, Sweden. Phone: +46 (0)346 124 25 E-mail: post@granberg-ab.se

User Manual issue date: 17.01.2023

granberggloves.com

EN USER MANUAL FOR DISPOSABLE GLOVES CATEGORY III and MEDICAL DEVICE

The User Manual should be used with product-specific information.

User Instructions should be read before using.

INTENDED USE

Powder-free examination and protective disposable nitrile gloves are intended for use in the medical field to protect patients and users from cross-contamination. These gloves are also intended to protect against certain chemicals, microorganisms where hand protection is needed.

WARNINGS AND PRECAUTIONS OF USE

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals and other factors influencing the performance such as temperature, abrasion, degradation etc. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested.

PRODUCT INSTRUCTION FOR USE

Before use, after donning, and during use inspect the gloves for any defect or imperfections and discontinue use immediately if signs of tearing, swelling or degradation, or any damage appear. Dry hand before donning. Ensure chemicals or residuals cannot enter through the cuff. Change glove after each patient. Always select the correct size glove for your hand.

DISPOSAL

Used gloves can be contaminated and must be disposed of under hospital policy and/or local regulation.

INGREDIENTS/HAZARDOUS COMPONENTS

Components used in glove manufacturing may cause allergic reactions in some users. If allergic reactions occur, seek medical advice immediately. Where relevant, a list of substances contained in the glove that are known to cause allergies, per listed in Annex G of EN ISO 21420:2020, shall be supplied on request.

STORAGE

Store in a cool and dry place in its original package. Opened boxes should be kept away from fluorescent and sunlight. Keep the gloves away from ozone, heating devices, and the source of the fire. Gloves are packed in a dispenser box suitable for transport. Keep the gloves in the box when not in use.

REPORTING OF INCIDENTS

In case of any serious incident occurred with the use of this device, please report it to the manufacturer and the competent Authority.

Further information can be obtained from the manufacturer, please contact Granberg AS.

EXPLANATION OF SYMBOLS AND PICTOGRAMS USED

Protective gloves against dangerous chemicals and microorganisms - Part 1: Terminology and performance requirements for chemical risks. EN ISO 374-1:2016+A1:2018. Definition of breakthrough time through the glove palm (1 µg/cm²/min). Type A > level 2 for 6 chemicals, Type B > level 2 for 3 chemicals, Type C > level 1 for 1 chemical (no code under pictogram).


Table with columns: ISO 374-1 Type A, B, C, A-B-C-D-E-F-G-H-I, A: Methanol, B: Acetone, C: Acetonitril, D: Dichloromethane, E: Carbon disulfide, F: Toluene, G: Diethylamine, H: Tetrahydrofuran, I: Ethyl acetate, J: n-Heptane, K: Sodium hydroxide 40%, L: Sulphuric acid 96%, M: Nitric acid 65%, N: Acetic acid 98%, O: Ammonium hydroxide 25%, P: Hydrogen peroxide 30%, S: Hydrofluoric acid 40%, T: Formaldehyde 37%

Additional information on chemical resistance obtainable from manufacturer.

Table with columns: Permeation Performance Level, Measured Breakthrough Time (minutes)

*Indicates that the glove falls below the minimum performance level as stated in EN ISO 374-1:2016+A1:2018 for the given individual hazard.

ISO 374-6:2016


Protection against bacteria, fungi and viruses

VIRUS


Keep away from sunlight

ISO 374-6:2016

Protection against bacteria and fungi, not tested against viruses

ISO 374-6:2016

Suitable for contact with foodstuffs. Note: not all gloves that are suitable for handling food may be suitable for all types of food. Check the Food Declaration of Compliance.

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