

POWDER-FREE BLUE NITRILE GLOVE



Product Presentation

Package of 500 Units

- Contents: 10 boxes of 50 units
- **Dimensions:** 340 x 260 x 270 mm

Box of 50 Units

- Contents: 50 units
- Dimensions: 260 x 125 x 65 mm

Labelling

- Name and address of the Manufacturer
- Product name in different languages
- Commercial reference, Batch and bar code Storage conditions
- AQL 1.5
- Expiry
- CE Marking

- Size and number of units
- Disposable
- Protection pictograms
- Laws and Standards of reference
- Uses, applications and warnings

General Characteristics

Description: Extra-strong and extra-long powder-free blue nitrile glove.

Thin and extra sensitive to the touch due to the finger area being textured for better grip in dry and wet conditions. Reinforced rolled cuff.

The surface of the glove is treated with chlorine to prevent the gloves from sticking together and for easier on and off.

The nitrile offers three times greater protection against the microorifices present in conventional latex gloves, making it the best choice for a latex-free glove.

Useful life: 5 years





POWDER-FREE BLUE NITRILE GLOVE

🔕 TECHNICAL DATA 🖞 🚦 💽

Classification:

Medical Device **Class I**; *Royal Decree 1591/2009, Regulation (EU) 2017/745.* PPE **Category III**; *Regulation (EU) 2016/425*

Sizes: Small, Medium, Large, Extra Large

Colour: Blue.

| Property Performance level/Result | | Standards applied | | | |
|--|---|--|---------------------------|--|--|
| | Medical devices Regulation (EU | | | | |
| Absence of holes | Compliant | EN 455-1 | | | |
| Dimensions | Compliant | EN 455-2 | | | |
| Breaking strength | Compliant (average=13 N) | | | | |
| Biological safety requirements. | Compliant | EN 455-3 | | | |
| Breaking tensile strength | 13 Mpa | | | | |
| Elongation | 500% | | | | |
| | ASTM Sta | indards | | | |
| NBR medical gloves | Compliant | ASTM D6319 - 10 | | | |
| Residual powder | Compliant | ASTM D6124 - 06 | | | |
| Permeability to cytostatics | Compliant ASTM D6978 - 05 | | | | |
| Personal Protective Equipment Regulation (EU) 2016/425 | | | | | |
| Dexterity | 5 | | EN 420:2003+A1:2009 | | |
| Resistance to permeati | on by micro-organisms | T | | | |
| Air leak test | Compliant EN 374-2:2016 | | | | |
| Water leak test | Compliant | EN 374-5:2016 | | | |
| | | | VIRUS | | |
| Resistance to permeation by chemical products | | | | | |
| (J) n-Heptane | Class 1/ Permeation time >10 min | | | | |
| (K) Sodium Hydroxide (40%) | Class 6/ Permeation time >480 min | КРТ | | | |
| Methanol in water (1.5%) | Class 6/ Permeation time >480 min | EN 374-1:2016 | | | |
| Phenol (0.1%) | Class 6/ Permeation time >480 min | EN 16523-1: 2015 | | | |
| Paseo de la Castellana no. 141 Plantas 18 y 19 28046 Madrid Tel: 902367762 | <u>www.santex.es</u> <u>calidad@santex.es</u> <u>santex@santex.es</u> | Written by: <i>Magda Ridorsa</i> Latest version 21/07/2022 | Page 2 of 5 | | |



GD24 POWDER-FREE BLUE NITRILE GLOVE 🛞 TECHNICAL DATA 🖞 🛊 💽

| Glutaraldehyde (50%) | Class 6/ Permeation time >480 | | | | |
|---|--|--|--|--|--|
| | min | | | | |
| (P) Hydrogen | Class 5 / Permeation time > | | | | |
| peroxide (30%) | 240 min | | | | |
| (T) Formaldehyde | Class 6/ Permeation time >480 | | | | |
| (37%) | min | | | | |
| Ethidium Bromide | Class 6/ Permeation time >480 | | | | |
| (5%) | min | | | | |
| Plastic materials intended to come into contact with food Regulation ((EU) 2016/1416) | | | | | |
| Plastic materials inte | nded to come into contact with fo | od Regulation ((EU) 2016/1416) | | | |
| Plastic materials inter Migration | nded to come into contact with fo | od Regulation ((EU) 2016/1416) EN 1186-7:2002 | | | |
| | nded to come into contact with fo Compliant | | | | |
| Migration | | EN 1186-7:2002 | | | |
| Migration - Acetic | Compliant | EN 1186-7:2002 | | | |
| Migration - Acetic - Ethanol 10% /20% | Compliant | EN 1186-7:2002 | | | |

Physical Properties

Composition: Nitrile Butadiene Rubber (NBR)

Characteristics:

- Ambidextrous

- AQL: 1.5
- Latex-Free



- Powder-Free
- Chlorinated Free from Proteins and chemical accelerants
- Textured Fingers
- Free from Thiuram, animal tissue and other biological substances

Sizes

| Size | Weight | Length | Palm | Thickness (mm) ±0.02 | | |
|------|---------|--------|-------------------|----------------------|------|-------------------|
| | (g)±0.3 | (mm) | width (mm) ± 5 | Finger | Palm | Cuff (25 cm long) |
| S | 8.4 | ≥285 | 85 | 0.17 | 0.11 | 0.08 |
| Μ | 8.8 | ≥285 | 95 | 0.17 | 0.11 | 0.08 |
| L | 9.2 | ≥285 | 106 | 0.17 | 0.11 | 0.08 |
| XL | 9.6 | ≥285 | 116 | 0.17 | 0.11 | 0.08 |

| Paseo de la Castellana no. 141 Plantas 18 y 19 28046 Madrid Tel: 902367762 | www.santex.es calidad@santex.es santex@santex.es | Written by: <i>Magda Ridorsa</i> Latest version 21/07/2022 | Page 3 of 5 |
|--|--|--|---------------------------|
|--|--|--|---------------------------|



POWDER-FREE BLUE NITRILE GLOVE

TECHNICAL DATA 🖞 👔

Logistics sheet

| REF - Size | EAN Inner box | Code Package | Kg Packages | Volume m ³ | Boxes/ Pallet | Assembly/ Pallet (Boxes x heights) |
|------------|---------------|-----------------|----------------|--------------------------|------------------|---|
| GD24BB- S | 8437001266708 | 8437001266746 | 5.10 | 0.023868 | 63 | 9x7 |
| GD24BC- M | 8437001266715 | 8437001266753 | 5.45 | 0.023868 | 63 | 9x7 |
| GD24BD- L | 8437001266722 | 8437001266760 | 5.60 | 0.023868 | 63 | 9x7 |
| GD24BE- XL | 8437001266739 | 8437001266777 | 6.00 | 0.023868 | 63 | 9x7 |

Uses and applications

In the **medical** field, gloves for doing **medical examinations**, orthodontics, clinical examinations, diagnostic and therapeutic procedures, for laboratory use, and in general for all activities which require gloves to create a protective barrier against infectious bodies, like in the fields of research and veterinary science.

They are suitable for all uses requiring high performance. They protect users against contamination from potentially infectious materials, such as viruses, bacteria, blood and infected body fluids and other contaminating materials. They offer low protection against **chemical** risks. They comply with requirements for confirming microbiological safety and low chemical risk (EN374-2 and EN374-4).

They are also used in the **food**, **electronics** and **cleaning** industries, because NBR does not contain latex or chemical accelerants, so they reduce skin irritation due to allergies and provide acceptable comfort and elasticity. In the food sector, these gloves comply with the requirements of Regulation 10/2011 regarding plastic materials intended to come into contact with food.

Storage Conditions

Store in a cool dry place. Avoid excess heat and protect from direct sunlight and fluorescent lighting.





POWDER-FREE BLUE NITRILE GLOVE



Directives and Standards of reference

- EN 374/1-2-4-5: Protective gloves against chemicals and micro-organisms.
- EN 420: Protective gloves, general requirements and testing methods.
- EN 1186/7:2002: Testing methods for overall migration into aqueous food simulants using a pouch.
- EN 455/1-2-3: Disposable medical protective gloves.
- ISO 13485:2003: Quality System for the manufacturing of Medical Devices.
- ASTM D 6124-06: Residual powder content.
- ASTM D 6319: Standard specification for nitrile examination gloves for medical use.
- ASTM D 6978-05: Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.
- Regulation 10/2011 on plastic materials and articles intended to come into contact with foo Text relevant for the purposes of the EEE.
- Royal Decree 866/2008, regarding plastic materials intended to come into contact with food which transposes directive 2002/72/EEC, repealed by Regulation 10/2011.
- Regulation (EU) 2017/745, which regulates Medical Devices.
- Regulation (EU) 2020/1245, which modifies and corrects Regulation (EU) 10/2011 on plastic materials and objects intended to come into contact with food.

Management System

Management System compliant with the ISO 13485 standards.



| Paseo de la Castellana no. 141 Plantas 18 y 19 28046 Madrid Tel: 902367762 | www.santex.es calidad@santex.es santex@santex.es | Written by: <i>Magda Ridorsa</i> Latest version 21/07/2022 | Page 5 of 5 |
|--|--|--|---------------------------|
|--|--|--|---------------------------|